

ClearSkin Hub

iPad Therapist Workflow Guide

A short step-by-step training guide for daily therapist use on iPad, from install to scheduling, consent, sim worksheet, and note generation.

Start here on iPad Safari: <https://dermathetherapies.vercel.app/>

All screenshots in this guide use demo patients created for training. No real PHI is shown.

ClearSkin Hub

Set a local PIN to protect patient notes on this device.

Add to Home Screen First

Do this first before setup, or you may need to enter everything again.

1. Tap the Safari Share button, the square with an upward arrow.
2. Tap View More, then scroll down and tap Add to Home Screen.
3. Tap Add.

Safari only. Add it to Home Screen before creating your PIN. [Maybe Later](#)

Default Therapist

Supervising Physician Name

Dermatology Office Name

Dermatology Office Logo



Set PIN

[Save Setup](#)

Figure: First-time iPad setup screen with Add to Home Screen guidance.

Table of Contents

Use this guide in order the first time, then keep it as a quick reference during daily use.

1. Add ClearSkin Hub to the iPad Home Screen	3
2. Learn the Main Navigation	4
3. Use Schedule to Plan Visits	5
4. Add a New Patient	6
5. Use the Patient Chart as the Main Record	7
6. Start with Consent / Path Intake	8
7. Review and Sign Consent	9
8. Complete Course Setup and Generate the Sim Worksheet	10
9. Access Consent and Sim Documents from the Course	11
10. Generate Sim / Consult or Treatment Notes	12
11. Enter Vitals, Photos, and Attachments on iPad	13
12. Restore Unsaved Edits After iPad Idle Time	14
13. Use Active Patients for Daily Treatment	15
14. Use the Consent / Sim Docs Tab for Document-Only Work	16
What Therapists Can Do In ClearSkin Hub	17

1. Add ClearSkin Hub to the iPad Home Screen

- In Safari on the iPad, go to <https://dermatherapies.vercel.app/> first.
- Open ClearSkin Hub in Safari only.
- Tap Share, then View More, then Add to Home Screen before doing first-time setup.
- Add it once so the app opens like a full-screen iPad app during daily use.
- On first setup, enter the default therapist, supervising physician, office name, and a 4-8 digit PIN.

After setup, the app shows a recovery code one time only. Write it down and store it offline.

CLEAR SKIN CANCER

ClearSkin Hub

Set a local PIN to protect patient notes on this device.

Add to Home Screen First

Do this first before setup, or you may need to enter everything again.

1. Tap the Safari Share button, the square with an upward arrow.
2. Tap View More, then scroll down and tap Add to Home Screen.
3. Tap Add.

Safari only. Add it to Home Screen before creating your PIN. [Maybe Later](#)

Default Therapist

Supervising Physician Name

Dermatology Office Name

Dermatology Office Logo

Set PIN

New PIN


Confirm PIN

[Save Setup](#)

The first setup screen includes the Home Screen install reminder.

2. Learn the Main Navigation

- Schedule opens the weekly calendar for consults, treatments, follow-ups, clinic hours, and closed time.
- Active Patients is the day-to-day treatment list for patients currently moving through care.
- Completed Patients keeps finished courses easy to reopen and review.
- Archive hides inactive patients from the daily list without deleting history.
- Consent / Sim Docs is for document-only workflows when you need consent and worksheet handling outside a full active course.



Main navigation after setup.

3. Use Schedule to Plan Visits

- Tap Schedule to view the Monday-Friday clinic week.
- Use Previous Week, Current Week, and Next Week to move through the calendar. The week field shows the first and last workday.
- Use Clinic Hours to set regular open/close times and add closed office time or holidays.
- Tap Add Appointment, or tap a calendar time slot, to add sim/consult, treatment, or follow-up visits.
- Manual appointments can be entered without using the patient chart side of the app.
- Linked course treatment schedules use projected or prescribed fractions, appointment length, and selected weekdays to create the full treatment schedule automatically.
- A time slot can hold up to two appointments. A third overlapping appointment is blocked so the calendar does not get overbooked.
- Tap an appointment to start the scheduled note, update the status, or edit the appointment. Drag an appointment when only the date or time needs to change.
- Print Schedule is available from the pending course so the patient can receive the calendar before the sim worksheet is completed.

When a pending sim/consult appointment is started from the calendar, it can move into the active patient workflow without retyping the intake details.

4. Add a New Patient

- Tap Add Patient from Active Patients.
- Enter first name, last name, MRN, DOB, and sex.
- Use the on-screen number pad for DOB entry on iPad.
- Tap Save Patient to open the chart.

ClearSkin Hub

Local-first treatment note workflow

Active Workflow

Add Patient

First Name

Last Name

MRN

DOB

Sex


Cancel

Active Patients

Patient creation modal on iPad.

5. Use the Patient Chart as the Main Record

- The patient chart keeps demographics, courses, notes, PDFs, and linked documents together.
- Tap Add Course to start a new clinical course.
- Tap Edit Patient when demographics need to be corrected.
- Archive Patient removes the chart from the active list but keeps the full history available.



Newly created patient chart before any course is added.

6. Start with Consent / Path Intake

- Use intake to capture what is already known before full treatment setup is finished.
- Enter lesion count, biopsy date, sim / consult date, treatment location, diagnosis, ICD10, and projected fractions.
- Tap Save Intake to create a pending course.
- This keeps the patient moving forward even if final setup details are not ready yet.

ClearSkin Hub

Local...

New Consent / Path Intake

Start the course with the pathology details we already know. Full treatment setup can be completed after the sim / consult.

Projected Fractions are an estimate for the Sim Worksheet and can be changed on the 1st treatment to the Final Prescribed Fractions.

Number of Lesions: 1 Lesion

Biopsy Date: 04/21/2026

Sim / Consult Date: 04/21/2026

Lesion

Treatment Lesion: Treatment location

Diagnosis: Select Diagnosis

ICD10: ICD10

Projected Fractions: Select Fractions

Cancel Save Intake

Carson, Maya

MRN 100245 · DOB 07/12/1968

Edit Patient Add Course Archive Patient

Consent / Path Intake starts the course with only the essentials.

7. Review and Sign Consent

- From the pending course, tap Review / Sign Consent.
- Read the form, capture the patient signature, then capture the witness signature.
- Once saved, the consent PDF stays linked to the course.
- If a signed paper form already exists, use Import or Replace instead of re-signing.

Review Consent Form

Carson, Maya
MRN 100245 · DOB 1968-07-12
1-lesion course

Sites
Lesion 1: Left nasal ala · Basal Cell Carcinoma

Permission Granted: I agree to undergo radiation therapy for the treatment of my skin cancer. I am aware of alternative medical and surgical treatments, as well as the risks and benefits of radiation therapy.

Acknowledgment Necessity to Complete Course of Treatment: I understand that the effectiveness and potential cure rate of my prescribed radiation treatment depends on my full participation in the treatment plan. This includes attending all scheduled appointments and receiving each required dose of radiation as prescribed by my provider. I acknowledge that missed or delayed treatments may reduce the likelihood of achieving the intended therapeutic outcome.

Immobilization and Photographs: The use of photographs for treatment purposes has been explained to me. I give permission, as part of my treatment, to have photographs taken as needed. Treatments may require the use of temporary immobilization. I give permission to use such devices as needed.

Risks & Complications Explained: I have been made aware of risks, such as but not limited to redness, dryness, peeling of skin, possible bleeding, risks of infection, scarring, radiation dermatitis, prolonged wound healing, nerve injury, inability to clear the tumor, possible permanent hypo/hyper pigmentation, telangiectasia, and possible permanent hair loss in immediate treatment area. I have been informed there is a remote possibility of secondary cancers developing in the treated area in the future.

No Guarantee: I acknowledge that no guarantee of complete success has been made to me concerning the result intended from the therapy.

Understanding This Form: I confirm that I have read and fully understand this form and that all blank spaces have been completed or crossed off prior to my signing. I have been given an opportunity to ask questions, and all my questions have been answered fully and to my satisfaction.

☐ I have read and understand this consent form and am ready to proceed to signatures.

Cancel Continue To Patient Signature

1-lesion course · pending Edit Intake Review / Sign Consent Import Signed Consent Complete Course Setup

Lesion 1

Consent review is step-based so signatures happen in order.

8. Complete Course Setup and Generate the Sim Worksheet

- After consent/path intake, tap Complete Course Setup.
- Enter the treatment-specific details such as cone, cutout, lesion size, depth, positioning, side, and any setup devices.
- Tap Save Course to activate the course.
- Saving the completed setup automatically generates the Sim Worksheet.

Complete Course Setup

Face Photo
Choose File No file chosen

Add if one is not already on file.

Lesion
Treatment Lesion: Left nasal ala
Diagnosis: Basal Cell Carcinoma
ICD10: C44.311

Cone Size
Select Cone ▼

Cutout Size
Open Cone ▼

Lesion Size (mm)
e.g. 10mm

Treatment Depth
3 mm ▼

Additional Treatment Devices

- ☐ Eye Shield
- ☐ Ear Shield
- ☐ Nasal Shield
- ☐ Gum Shield
- ☐ Hand Ring
- ☐ Knee Wedge
- ☐ Pillow
- ☐ Flash Shield
- ☐ Lip Shield
- ☐ Special Set-up

Positioning

- ☐ Supine
- ☐ Prone
- ☐ Vac-Lok
- ☐ Sitting in Chair
- ☐ Hunched Over Tx Couch

Lesion Side

- ☐ Left
- ☐ Right
- ☐ Medial


Cancel Remove Course Save Course

1-lesion course · pending Intake Consent Form Re-sign Consent Replace Consent Complete Course Setup

The full setup screen finishes the active treatment course.

9. Access Consent and Sim Documents from the Course

- On the patient chart, tap Documents on the course row.
- Open Consent Form launches the linked signed consent PDF.
- Open Sim Worksheet launches the worksheet PDF created during setup.
- Re-sign, replace, or regenerate from the same window if details change later.
- Generated consent and sim worksheet PDFs open in a preview first so you can confirm the correct form before using the iPad share button to save it.
- When saved from the preview, the file name should use the patient name, such as Patient Name - Consent or Patient Name - Sim Worksheet.



Each course keeps its consent and sim worksheet together.

10. Generate Sim / Consult or Treatment Notes

- Tap Start Today's Note from the active course.
- The note screen opens with Visit Details, Note Preview, Save Draft, and Finalize.
- Set Visit Type to Sim / Consult or treatment as needed.
- Sim / Consult notes include Tx site name, ICD10, Total Fractions, Daily dose, Prescribed dose, Tx depth, Cone size, Cutout flex shield size, and Additional Tx devices.
- Daily dose and prescribed dose are filled from the projected or prescribed fraction logic and display in cGy.
- Save Draft keeps the note editable. Finalize creates the current PDF output for the visit.

Carson, Maya

Left nasal ala · 1-lesion course

Visit Type: Sim / Consult

Visit Details

Visit Date: 04/21/2026

Visit Type: Sim / Consult

Therapist: Jordan Ellis RT(T)

Projected Fractions: 10

Additional Notes

Biopsy Date: 04/21/2026

Treatment Start Date: MM/DD/YYYY

☒ Include Exam Vitals

Exam Vitals

Blood Pressure: e.g. 120/80

Heart Rate: e.g. 72 BPM

O2 Saturation: e.g. 98%

Weight: e.g. 165 lbs

☐ Ultrasound Performed ☐ Add MIPS

Visit editor with save actions across the top.

11. Enter Vitals, Photos, and Attachments on iPad

- Use the docked number pad for dates, blood pressure, lesion size, heart rate, O2 saturation, weight, and PIN entry.
- Units like mm, BPM, %, lbs, and mmHg are added automatically where supported.
- Attach daily treatment photos and extra files directly from the visit.
- Linked course documents are available lower in the visit so everything stays connected.

The screenshot displays the iPad interface for a patient named Carson, Maya. At the top, the patient's name is shown, followed by the treatment description 'Left nasal ala · 1-lesion course'. Action buttons for 'Visit Details', 'Note Preview', 'Save Draft', and 'Finalize' are present. The 'Visit Type' is set to 'Sim / Consult'. The 'Visit Details' section includes fields for 'Visit Date' (04/21/2026), 'Visit Type' (Sim / Consult), and 'Therapist' (Jordan Ellis RT(T)). There is also a 'Projected Fractions' dropdown set to 10 and an 'Additional Notes' text area. Below this, 'Biopsy Date' (04/21/2026) and 'Treatment Start Date' (MM/DD/YYYY) are entered. A checkbox for 'Include Exam Vitals' is checked. The 'Exam Vitals' section contains fields for 'Blood Pressure' (e.g. 120/80), 'Heart Rate' (e.g. 72 BPM), 'O2 Saturation' (e.g. 98%), and 'Weight' (e.g. 165 lbs). At the bottom, there are checkboxes for 'Ultrasound Performed' and 'Add MIPS'. A 'Lesion 1' section lists treatment details: 'Treatment Lesion: Left nasal ala', 'Diagnosis: Basal Cell Carcinoma', 'ICD10: C44.311', 'Treatment Depth: 3 mm', 'Daily Dose: -', and 'Total Dose: -'.

Carson, Maya

Left nasal ala · 1-lesion course

Visit Type: Sim / Consult

Visit Details

Visit Date: 04/21/2026 Visit Type: Sim / Consult Therapist: Jordan Ellis RT(T)

Projected Fractions: 10

Additional Notes

Biopsy Date: 04/21/2026 Treatment Start Date: MM/DD/YYYY

☒ Include Exam Vitals

Exam Vitals

Blood Pressure: e.g. 120/80 Heart Rate: e.g. 72 BPM O2 Saturation: e.g. 98%

Weight: e.g. 165 lbs

☐ Ultrasound Performed ☐ Add MIPS

Lesion 1

Treatment Lesion: Left nasal ala

Diagnosis: Basal Cell Carcinoma

ICD10: C44.311

Treatment Depth: 3 mm

Daily Dose: -

Total Dose: -

The iPad note workflow is optimized for fast number entry and attachment handling.

12. Restore Unsaved Edits After iPad Idle Time

- If the iPad has been sitting in the background for a long time, Safari may reload or resume the app in the middle of a workflow.
- ClearSkin Hub keeps recovery drafts for major forms, including new/edit patient, course/intake/setup, consent signing, document-only records, and sim worksheet setup.
- If unsaved work is found, the app shows an Unsaved Edits Found prompt.
- Tap Restore Edits to bring the typed information back into the form.
- Tap Discard only when you are sure you no longer need those unsaved edits.
- A successful Save, Finalize, or an intentional Cancel clears the recovery draft so it does not keep appearing.

Recovery drafts are for typed form information and signatures. Large uploaded files or photos are not stored in recovery drafts to avoid iPad storage issues.

13. Use Active Patients for Daily Treatment

- Once the course is active, it shows up in Active Patients with quick actions.
- Start Today's Note opens the next note for the course.
- Resume Last Note appears when a draft already exists.
- When the course is finished, move it to Completed Patients without losing history.

The screenshot displays the patient profile for **Carson, Maya** with MRN 100245 and DOB 07/12/1968. At the top right are buttons for 'Edit Patient', 'Add Course', and 'Archive Patient'. Below the patient name, the course 'Left nasal ala' is shown with 'Sim / Consult' and '1-lesion course · active' status. Action buttons for 'Documents', 'Edit Course', and 'Start Today's Note' are present. A section titled 'Lesion 1' contains the text 'Left nasal ala', 'C44.311', and 'Dose set on 1st treatment'. A blue notification banner at the bottom right states 'Course setup completed and sim worksheet generated.'

Active course row after setup is complete.

14. Use the Consent / Sim Docs Tab for Document-Only Work

- Use this tab when you need consent and sim worksheet handling without building a full active treatment course first.
- Tap Add Patient Info to create a document-only record.
- This is useful for early paperwork, standalone documentation, or worksheet prep before the full course is ready.
- For active treatment patients, the patient chart remains the best place to manage documents.



Document-only workflow starts from the Consent / Sim Docs tab.

What Therapists Can Do In ClearSkin Hub

- Create patients, run one-lesion or two-lesion courses, and keep finished history searchable.
- Build clinic schedules for sim/consult, treatment, and follow-up visits, either by linked course or manual appointment.
- Capture consent, sim worksheet setup, treatment vitals, note text, daily photos, and extra attachments in one chart.
- Reopen an old note, correct the structured fields or text, and regenerate the PDF without creating a duplicate visit.
- Recover unsaved form work after the iPad has been idle, reloaded, or resumed from the background.
- Use iPad-friendly number pads for PINs, dates, blood pressure, lesion size, heart rate, O2 sat, and weight.
- Keep all data local to the device or local app storage. The app is designed for local-first use, not cloud workflow.

Best daily path

Open **Schedule** for the clinic day, or open **Active Patients** -> open patient -> use **Documents** when consent/worksheet needs attention -> tap **Start Today's Note** -> save draft or generate PDF.